

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 050

**Graduated cylinder, 500:5 ml, A, conf. cert.**




**Batch:** 30-21

**Serial:**

**Measuring Nr.** 48091

**Batch size:** 795

Volume:	100,000 ml	300,000 ml	500,000 ml
<b>Tolerance:</b>	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
<b>Average value:</b>	99,9988 ml	300,1681 ml	500,3816 ml
<b>Deviation:</b>	0,3854 ml	0,3301 ml	0,2789 ml
<b>Precision (CV):</b>	0,3854 %	0,1100 %	0,0558 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	29.07.2021
<b>Balance:</b>	Sartorius - LA 1200 S

