

Certificate of performance



Article: 3 715 050

Graduated cylinder, 500:5 ml, A, conf. cert.




Batch: 44-22

Serial:

Measuring Nr. 60271

Batch size: 280

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
Average value:	100,2913 ml	300,0905 ml	499,6846 ml
Deviation:	0,4601 ml	0,2860 ml	0,3273 ml
Precision (CV):	0,4601 %	0,0953 %	0,0655 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	11.11.2022
Balance:	Sartorius - LA 1200 S

