

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050

Graduated cylinder, 500:5 ml, A, conf. cert.




Batch: 36-23

Serial:

Measuring Nr. 64415

Batch size: 1200

| Volume: | 100,000 ml | 300,000 ml | 500,000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 2,5000 ml | ± 2,5000 ml | ± 2,5000 ml |
| Average value: | 100,0414 ml | 300,2567 ml | 500,1363 ml |
| Deviation: | 0,2893 ml | 0,2220 ml | 0,2370 ml |
| Precision (CV): | 0,2893 % | 0,0740 % | 0,0474 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|--------------------------|
| Operator: | 7 |
| Date: | 15.09.2023 |
| Balance: | Sartorius - LA 1200 S |

