

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 100 BL USP

Grad. cylinder, 1000:10 ml, blue, USP, A, conf. cert.




Batch: 08-19

Serial: -

Measuring Nr.: 27981

Batch size: 30

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 3,0000 ml	± 3,0000 ml	± 3,0000 ml
Average value:	199,1245 ml	599,4358 ml	1000,0299 ml
Deviation:	0,7519 ml	0,5268 ml	0,3508 ml
Precision (CV):	0,3760 %	0,0878 %	0,0351 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	25.02.2019
Balance:	Sartorius - LC 6200 S

