

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 100 BL

**Grad. cylinder, 1000:10 ml, blue grad., A, conf. cert.**




**Batch:** 34-18

**Serial:**

**Measuring Nr.:** 24143

**Batch size:** 280

| Volume:                | 200,000 ml  | 600,000 ml  | 1000,000 ml  |
|------------------------|-------------|-------------|--------------|
| <b>Tolerance:</b>      | ± 5,0000 ml | ± 5,0000 ml | ± 5,0000 ml  |
| <b>Average value:</b>  | 200,5601 ml | 600,3149 ml | 1000,8298 ml |
| <b>Deviation:</b>      | 0,8184 ml   | 0,6543 ml   | 0,8785 ml    |
| <b>Precision (CV):</b> | 0,4092 %    | 0,1091 %    | 0,0879 %     |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                          |
|------------------|--------------------------|
| <b>Operator:</b> | 4                        |
| <b>Date:</b>     | 17.09.2018               |
| <b>Balance:</b>  | Sartorius - LA<br>1200 S |

