

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 100 BL

Grad. cylinder, 1000:10 ml, blue grad., A, conf. cert.

Batch: 06-23




Serial:

Measuring Nr. 61647

Batch size:

237

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	200,0811 ml	600,0401 ml	999,5704 ml
Deviation:	0,7602 ml	0,6609 ml	0,7883 ml
Precision (CV):	0,3801 %	0,1101 %	0,0788 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

09.02.2023

Balance:

Sartorius - LC
6200 S

