

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 100 USP

**Graduated cylinder, 1000 ml, USP, A, conf. cert.**




**Batch:** 24-21

**Serial:**

**Measuring Nr.:** 47581

**Batch size:** 112

| Volume:                | 200,000 ml  | 600,000 ml  | 1000,000 ml  |
|------------------------|-------------|-------------|--------------|
| <b>Tolerance:</b>      | ± 3,0000 ml | ± 3,0000 ml | ± 3,0000 ml  |
| <b>Average value:</b>  | 200,5206 ml | 600,0791 ml | 1000,3002 ml |
| <b>Deviation:</b>      | 0,5582 ml   | 0,8339 ml   | 0,3149 ml    |
| <b>Precision (CV):</b> | 0,2791 %    | 0,1390 %    | 0,0315 %     |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                          |
|------------------|--------------------------|
| <b>Operator:</b> | 4                        |
| <b>Date:</b>     | 23.06.2021               |
| <b>Balance:</b>  | Sartorius - LC<br>6200 S |

