

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 100 USP

**Graduated cylinder, 1000 ml, USP, A, conf. cert.**




**Batch:** 39-21

**Serial:**

**Measuring Nr.:** 49342

**Batch size:** 178

| Volume:                | 200,000 ml  | 600,000 ml  | 1000,000 ml |
|------------------------|-------------|-------------|-------------|
| <b>Tolerance:</b>      | ± 3,0000 ml | ± 3,0000 ml | ± 3,0000 ml |
| <b>Average value:</b>  | 199,9645 ml | 600,0574 ml | 999,9094 ml |
| <b>Deviation:</b>      | 0,7802 ml   | 0,6715 ml   | 0,6340 ml   |
| <b>Precision (CV):</b> | 0,3901 %    | 0,1119 %    | 0,0634 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                          |
|------------------|--------------------------|
| <b>Operator:</b> | 4                        |
| <b>Date:</b>     | 28.09.2021               |
| <b>Balance:</b>  | Sartorius - LC<br>6200 S |

