

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 100 USP

Graduated cylinder, 1000 ml, USP, A, conf. cert.




Batch: 26-23

Serial:

Measuring Nr.: 63576

Batch size: 99

| Volume: | 200,000 ml | 600,000 ml | 1000,000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 3,0000 ml | ± 3,0000 ml | ± 3,0000 ml |
| Average value: | 200,6798 ml | 600,2285 ml | 999,9732 ml |
| Deviation: | 0,4090 ml | 0,7000 ml | 0,4046 ml |
| Precision (CV): | 0,2045 % | 0,1167 % | 0,0405 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|-------------------------------------|
| Operator: | 4 |
| Date: | 29.06.2023 |
| Balance: | Sartorius - MSE14202S- OCE-D0 |

