

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 100 USP

**Graduated cylinder, 1000 ml, USP, A, conf. cert.**




**Batch:** 09-24

**Serial:** -

**Measuring Nr.:** 66213

**Batch size:** 197

Volume:	200,000 ml	600,000 ml	1000,000 ml
<b>Tolerance:</b>	± 3,0000 ml	± 3,0000 ml	± 3,0000 ml
<b>Average value:</b>	200,2511 ml	600,4330 ml	1000,3969 ml
<b>Deviation:</b>	0,4396 ml	0,5279 ml	0,4303 ml
<b>Precision (CV):</b>	0,2198 %	0,0880 %	0,0430 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

04.03.2024

**Balance:**

Sartorius -  
MSE14202S-  
OCE-D0

