

Certificate of performance



Article: 3 715 100

Graduated cylinder, 1000:10 ml, A, conf. cert.




Batch: 38-21

Serial:

Measuring Nr.: 49406

Batch size: 415

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	199,7307 ml	599,9236 ml	999,8204 ml
Deviation:	0,9135 ml	0,6069 ml	0,6236 ml
Precision (CV):	0,4567 %	0,1011 %	0,0624 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	30.09.2021
Balance:	Sartorius - LC 6200 S

