

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 100

Graduated cylinder, 1000:10 ml, A, conf. cert.




Batch: 06-23

Serial:

Measuring Nr.: 61694

Batch size: 384

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	200,8267 ml	599,8395 ml	999,5262 ml
Deviation:	0,8633 ml	0,6695 ml	0,6010 ml
Precision (CV):	0,4317 %	0,1116 %	0,0601 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	14.02.2023
Balance:	Sartorius - LC 6200 S

