

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 200 BL

Grad. cylinder, 2000:20 ml, blue grad., A, conf. cert.




Batch: 45-21

Serial:

Measuring Nr. 50389

Batch size: 46

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,000 ml	± 10,000 ml	± 10,000 ml
Average value:	401,6974 ml	1198,9190 ml	1996,9488 ml
Deviation:	0,9124 ml	1,0309 ml	0,8373 ml
Precision (CV):	0,2281 %	0,0859 %	0,0419 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	11.11.2021
Balance:	Sartorius - LC 6200 S

