

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 200

**Grad. cylinders, 2000:20 ml, A, conf. cert.**




**Batch:** 36-17

**Serial:**

**Measuring Nr.** 13972

**Batch size:** 106

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,000 ml	± 10,000 ml	± 10,000 ml
Average value:	400,8796 ml	1200,7315 ml	1999,9310 ml
Deviation:	2,1229 ml	2,0680 ml	2,2008 ml
Precision (CV):	0,5307 %	0,1723 %	0,1100 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	07.09.2017
Balance:	Sartorius - LC 6200 S

