

Certificate of performance



Article: 3 715 200

Graduated cylinder, 2000:20 ml, A, conf. cert.




Batch: 03-18

Serial:

Measuring Nr. 17353

Batch size: 86

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,000 ml	± 10,000 ml	± 10,000 ml
Average value:	402,6619 ml	1201,0387 ml	1999,6708 ml
Deviation:	2,9257 ml	1,6456 ml	1,8016 ml
Precision (CV):	0,7314 %	0,1371 %	0,0901 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	22.01.2018
Balance:	Sartorius - LC 6200 S

