

Certificate of performance



Article: 3 715 200

Graduated cylinder, 2000:20 ml, A, conf. cert.




Batch: 31-18

Serial: -

Measuring Nr.: 23525

Batch size: 129

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,0000 ml	± 10,0000 ml	± 10,0000 ml
Average value:	399,2568 ml	1200,3658 ml	2000,5102 ml
Deviation:	1,5946 ml	1,5562 ml	2,2376 ml
Precision (CV):	0,3986 %	0,1297 %	0,1119 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	22.08.2018
Balance:	Sartorius - LC 6200 S

