

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 200

Graduated cylinder, 2000:20 ml, A, conf. cert.




Batch: 30-22

Serial:

Measuring Nr.: 57701

Batch size: 150

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,000 ml	± 10,000 ml	± 10,000 ml
Average value:	398,2245 ml	1198,4779 ml	1999,7913 ml
Deviation:	1,4382 ml	1,2870 ml	1,7573 ml
Precision (CV):	0,3595 %	0,1072 %	0,0879 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	28.07.2022
Balance:	Sartorius - LC 6200 S

