

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 003 BL

Grad. cylinder, 25:0,5 ml, blue grad., NS 14, conf. cert.




Batch: 24-19

Serial:

Measuring Nr.: 31040

Batch size: 78

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	4,9663 ml	14,9915 ml	24,9904 ml
Deviation:	0,0466 ml	0,0097 ml	0,0353 ml
Precision (CV):	0,9328 %	0,0644 %	0,1413 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	27.06.2019
Balance:	Sartorius - LE244S-OCE

