

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 754 003 BL

**Grad. cylinder, 25:0,5 ml, blue grad., NS 14, conf. cert.**




**Batch:** 08-24

**Serial:**

**Measuring Nr.** 66145

**Batch size:** 84

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	5,0134 ml	15,0192 ml	24,9953 ml
Deviation:	0,0177 ml	0,0250 ml	0,0394 ml
Precision (CV):	0,3534 %	0,1670 %	0,1577 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	21.02.2024
Balance:	Sartorius - LE244S-OCE

