

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 754 003

**Grad. cylinder, 25:0,5 ml, NS 14, conf.cert.**




**Batch:** 45-20

**Serial:**

**Measuring Nr.** 42273

**Batch size:** 150

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	4,9830 ml	15,0046 ml	25,0128 ml
Deviation:	0,0331 ml	0,0260 ml	0,0124 ml
Precision (CV):	0,6613 %	0,1736 %	0,0495 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	09.11.2020
Balance:	Sartorius - LE244S-OCE

