

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 003

Grad. cylinder, 25:0,5 ml, NS 14, conf.cert.




Batch: 08-24

Serial:

Measuring Nr. 66146

Batch size: 100

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	5,0081 ml	14,9726 ml	24,9532 ml
Deviation:	0,0138 ml	0,0107 ml	0,0221 ml
Precision (CV):	0,2767 %	0,0710 %	0,0886 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	22.02.2024
Balance:	Sartorius - LE244S-OCE

