

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 754 005

**Grad. cylinder, 50:1 ml, NS 19, A, conf. cert.**




**Batch:** 09-19

**Serial:**

**Measuring Nr.** 28175

**Batch size:** 150

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	10,0042 ml	30,0018 ml	50,0029 ml
Deviation:	0,1542 ml	0,0859 ml	0,0759 ml
Precision (CV):	1,5422 %	0,2863 %	0,1518 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	07.03.2019
Balance:	Sartorius - MSE224S-0CE -DU

