

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 010 BL

Grad. cylinders, 100:1 ml, NS 24, blue, A, conf. cert.




Batch: 10-19

Serial:

Measuring Nr. 28240

Batch size: 88

Volume:	20,000 ml	60,000 ml	100,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	20,0810 ml	60,0404 ml	100,0292 ml
Deviation:	0,1207 ml	0,0778 ml	0,0541 ml
Precision (CV):	0,6036 %	0,1297 %	0,0541 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	08.03.2019
Balance:	Sartorius - LA 1200 S

