

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 754 010 BL

**Grad. cylinders, 100:1 ml, NS 24, blue, A, conf. cert.**




**Batch:** 40-19

**Serial:**

**Measuring Nr.** 33485

**Batch size:** 90

Volume:	20,000 ml	60,000 ml	100,000 ml
<b>Tolerance:</b>	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
<b>Average value:</b>	19,9262 ml	59,9213 ml	100,0211 ml
<b>Deviation:</b>	0,1200 ml	0,1050 ml	0,1199 ml
<b>Precision (CV):</b>	0,6000 %	0,1751 %	0,1199 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	01.10.2019
<b>Balance:</b>	Sartorius - LE244S-OCE

