

Certificate of performance



Article: 3 754 010

Grad. cylinders, 100:1 ml, NS 24, A, conf. cert.




Batch: 13-21

Serial:

Measuring Nr. 45592

Batch size: 300

Volume:	20,000 ml	60,000 ml	100,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	19,9270 ml	59,9708 ml	100,0593 ml
Deviation:	0,0842 ml	0,0687 ml	0,0781 ml
Precision (CV):	0,4209 %	0,1144 %	0,0781 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	13.04.2021
Balance:	Sartorius - MSE224S-0CE -DU

