

Certificate of performance

Article: 3 754 025 BL

Grad. cylinder, 250:2 ml, NS 29, blue grad., A, conf. cert




Batch: 43-19

Serial:

Measuring Nr.: 34021

Batch size: 30

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
Tolerance:	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
Average value:	49,8891 ml	149,9495 ml	249,9394 ml
Deviation:	0,2298 ml	0,1210 ml	0,1786 ml
Precision (CV):	0,4595 %	0,0806 %	0,0715 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	31.10.2019
Balance:	Sartorius - LA 1200 S