

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 025 BL

Grad. cylinder, 250:2 ml, NS 29, blue grad., A, conf. cert




Batch: 49-23

Serial:

Measuring Nr. 65442

Batch size: 54

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
Tolerance:	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
Average value:	49,9308 ml	150,3038 ml	250,3987 ml
Deviation:	0,2735 ml	0,2167 ml	0,0601 ml
Precision (CV):	0,5470 %	0,1445 %	0,0240 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	07.12.2023
Balance:	Sartorius - LA 1200 S

DE-M