

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 754 025

**Graduated cylinder, 250:2 ml, NS 29, A, KB**




**Batch:** 49-23

**Serial:**

**Measuring Nr.:** 65477

**Batch size:** 52

Volume:	50,000 ml	150,000 ml	250,000 ml
<b>Tolerance:</b>	± 1,000 ml	± 1,000 ml	± 1,000 ml
<b>Average value:</b>	50,3192 ml	150,2573 ml	250,1417 ml
<b>Deviation:</b>	0,1100 ml	0,0893 ml	0,1141 ml
<b>Precision (CV):</b>	0,2200 %	0,0596 %	0,0456 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	12.12.2023
<b>Balance:</b>	Sartorius - LA 1200 S

DE-M