

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 050

Graduated cylinders, 500:5 ml, NS 34, A, conf. cert.




Batch: 30-19

Serial:

Measuring Nr.: 32021

Batch size: 40

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
Average value:	100,7607 ml	300,7626 ml	500,4253 ml
Deviation:	0,3029 ml	0,2532 ml	0,3020 ml
Precision (CV):	0,3029 %	0,0844 %	0,0604 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	07.08.2019
Balance:	Sartorius - LE244S-OCE

