

Certificate of performance



Article: 3 754 050

Graduated cylinders, 500:5 ml, NS 34, A, conf. cert.




Batch: 10-23

Serial:

Measuring Nr. 61946

Batch size: 104

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
Average value:	100,3369 ml	299,7808 ml	499,7294 ml
Deviation:	0,3367 ml	0,2431 ml	0,5071 ml
Precision (CV):	0,3367 %	0,0810 %	0,1014 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	14.03.2023
Balance:	Sartorius - LE244S-OCE

