

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 100 BL

Grad. cylinder, 1000 ml, NS 45, blue grad., A, conf. cert.




Batch: 49-17

Serial:

Measuring Nr. 16505

Batch size: 17

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	200,3940 ml	600,1714 ml	999,9079 ml
Deviation:	0,6470 ml	0,6917 ml	0,5608 ml
Precision (CV):	0,3235 %	0,1153 %	0,0561 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	15.12.2017
Balance:	Sartorius - LC 6200 S

