

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 100 BL

Grad. cylinder, 1000 ml, NS 45, blue grad., A, conf. cert.




Batch: 26-20

Serial:

Measuring Nr. 38877

Batch size: 30

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	199,6700 ml	599,6697 ml	1000,3190 ml
Deviation:	0,6136 ml	0,4380 ml	0,9832 ml
Precision (CV):	0,3068 %	0,0730 %	0,0983 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	30.06.2020
Balance:	Sartorius - LC 6200 S

