

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 754 100

**Grad. cylinder, 1000 ml, NS 45, A, conf. cert.**




**Batch:** 26-20

**Serial:**

**Measuring Nr.** 38973

**Batch size:** 44

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	200,4551 ml	600,2469 ml	1000,0047 ml
Deviation:	0,8126 ml	0,5878 ml	0,3524 ml
Precision (CV):	0,4063 %	0,0980 %	0,0352 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	02.07.2020
Balance:	Sartorius - LC 6200 S

