

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 754 200

**Grad. cylinder, 2000 ml, NS 45, A, conf. cert.**




**Batch:** 51-19

**Serial:**

**Measuring Nr.:** 35303

**Batch size:** 35

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,0000 ml	± 10,0000 ml	± 10,0000 ml
Average value:	400,2733 ml	1199,5115 ml	1998,7617 ml
Deviation:	2,4924 ml	2,4867 ml	1,8133 ml
Precision (CV):	0,6231 %	0,2072 %	0,0907 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	07.01.2020
Balance:	Sartorius - LC 6200 S

