

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 110

Grad. pipettes, 1:0,1 ml, AS, conf. cert.




Batch: 04-18

Serial:

Measuring Nr. 17703

Batch size: 1200

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
Average value:	0,9979 ml	0,5980 ml	0,1997 ml
Deviation:	0,0016 ml	0,0012 ml	0,0011 ml
Precision (CV):	0,1612 %	0,1963 %	0,5720 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	31.01.2018
Balance:	Sartorius - LE244S-OCE

DE-M