

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 642 001

Volumetric flask, 10 ml, GL 14, A, conf. cert.


Batch: 37-20

Serial: -

Measuring Nr. 40297

Batch size: 100

| | |
|------------------------|------------------|
| Volume: | 10,000 ml |
| Tolerance: | $\pm 0,0400$ ml |
| Average value: | 10,0040 ml |
| Deviation: | 0,0058 ml |
| Precision (CV): | 0,0583 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

21.09.2020

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M