

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 642 050

Volumetric flasks, 500 ml, GL 25, A, conf. cert.


Batch: 14-21

Serial:

Measuring Nr. 45479

Batch size: 20

Volume:	500,000 ml
Tolerance:	$\pm 0,2500$ ml
Average value:	500,0927 ml
Deviation:	0,0088 ml
Precision (CV):	0,0018 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	07.04.2021
Balance:	Sartorius - MSE14202S- OCE-D0

DE-M