

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 00214 BL

Volumetric flasks, 20 ml, NS 14, blue grad., A, conf. cer


Batch: 37-17

Serial:

Measuring Nr. 14316

Batch size: 500

Volume:	20,000 ml
Tolerance:	$\pm 0,0600$ ml
Average value:	19,9969 ml
Deviation:	0,0090 ml
Precision (CV):	0,0452 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	26.09.2017
Balance:	Sartorius - LE244S-OCE

