

# Certificate of performance

**Article:** 3 662 000 BL

**Vol. flasks, 5 ml, NS 7, blue grad., A, conf. cert.**


**Batch:** 37-18

**Serial:**

**Measuring Nr.** 24144

**Batch size:** 1168

<b>Volume:</b>	<b>5,0000 ml</b>
<b>Tolerance:</b>	$\pm 0,0250$ ml
<b>Average value:</b>	5,0005 ml
<b>Deviation:</b>	0,0051 ml
<b>Precision (CV):</b>	0,1018 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

9

**Date:**

17.09.2018

**Balance:**

Sartorius -  
LE244S-OCE