

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 00114 BL

Vol. flask, 10 ml, NS 14, blue grad., A, conf. cert.

Batch: 28-20


Serial:

Measuring Nr. 39091

Batch size:

315

Volume:	10,000 ml
Tolerance:	$\pm 0,0600$ ml
Average value:	10,0118 ml
Deviation:	0,0093 ml
Precision (CV):	0,0930 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	13.07.2020
Balance:	Sartorius - LE244S-OCE

DE-M