

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 662 002 USP

**Volumetric flasks, 20 ml, NS 10, USP, A, conf. cert.**


**Batch:** 46-23

**Serial:**

**Measuring Nr.** 65275

**Batch size:** 279

<b>Volume:</b>	20,000 ml
<b>Tolerance:</b>	± 0,0300 ml
<b>Average value:</b>	20,0059 ml
<b>Deviation:</b>	0,0030 ml
<b>Precision (CV):</b>	0,0152 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 7

**Date:** 24.11.2023

**Balance:** Sartorius -  
LE244S-OCE

DE-M