

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 002 W

Vol. flasks, 20 ml, NS 10, white grad., A, conf. cert.


Batch: 25-22

Serial:

Measuring Nr. 56964

Batch size: 200

| | |
|------------------------|-----------------|
| Volume: | 20,0000 ml |
| Tolerance: | $\pm 0,0400$ ml |
| Average value: | 19,9964 ml |
| Deviation: | 0,0045 ml |
| Precision (CV): | 0,0227 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 01.07.2022 |
| Balance: | Sartorius - LE244S-OCE |

DE-M