

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 003 BL

Vol. flasks, 25 ml, NS 10, blue grad., A, conf. cert.

Batch: 23-20


Serial:

Measuring Nr. 38794

Batch size:

3189

| | |
|------------------------|-----------------|
| Volume: | 25,0000 ml |
| Tolerance: | $\pm 0,0400$ ml |
| Average value: | 25,0027 ml |
| Deviation: | 0,0057 ml |
| Precision (CV): | 0,0229 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

18.06.2020

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M