

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 662 003

**Volumetric flasks, 25 ml, NS 10, A, conf. cert.**


**Batch:** 38-17

**Serial:**

**Measuring Nr.** 14302

**Batch size:** 1200

|                        |                   |
|------------------------|-------------------|
| <b>Volume:</b>         | <b>25,0000 ml</b> |
| <b>Tolerance:</b>      | $\pm 0,0400$ ml   |
| <b>Average value:</b>  | 24,9979 ml        |
| <b>Deviation:</b>      | 0,0065 ml         |
| <b>Precision (CV):</b> | 0,0259 %          |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                           |
|------------------|---------------------------|
| <b>Operator:</b> | 9                         |
| <b>Date:</b>     | 25.09.2017                |
| <b>Balance:</b>  | Sartorius -<br>LE244S-OCE |

DE-M