

Certificate of performance

Article: 3 662 005 BL

Vol. flasks, 50 ml, NS 12, blue grad., A, conf. cert.


Batch: 50-18

Serial:

Measuring Nr. 26617

Batch size: 2399

Volume:	50,000 ml
Tolerance:	$\pm 0,0600$ ml
Average value:	49,9953 ml
Deviation:	0,0096 ml
Precision (CV):	0,0192 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

20.12.2018

Balance:

Sartorius -
LE244S-OCE