

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 662 005

**Vol. flasks, 50 ml, NS 12, A, conf. cert**

**Batch:** 49-17


**Serial:**

**Measuring Nr.** 17661

**Batch size:**

1

|                        |                    |
|------------------------|--------------------|
| <b>Volume:</b>         | <b>50,000 ml</b>   |
| <b>Tolerance:</b>      | <b>± 0,0600 ml</b> |
| <b>Measured value:</b> | <b>50,0133 ml</b>  |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                           |
|------------------|---------------------------|
| <b>Operator:</b> | 9                         |
| <b>Date:</b>     | 30.01.2018                |
| <b>Balance:</b>  | Sartorius -<br>LE244S-OCE |

DE-M