

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 006 BL

Vol. flask, 50 ml, NS 14, blue grad., A, conf. cert.

Batch: 20-20


Serial:

Measuring Nr. 38351

Batch size:

1197

Volume:	50,0000 ml
Tolerance:	± 0,1000 ml
Average value:	49,9901 ml
Deviation:	0,0100 ml
Precision (CV):	0,0199 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

02.06.2020

Balance:

Sartorius -
LE244S-OCE

DE-M