

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 010 BL

Vol. flasks, 100 ml, NS 12, blue grad., A, conf. cert.

Batch: 08-23

Serial:

Measuring Nr. 61856

Batch size: 4985


Volume: 100,000 ml

Tolerance: $\pm 0,1000$ ml

Average value: 99,9982 ml

Deviation: 0,0154 ml

Precision (CV): 0,0154 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 7

Date: 02.03.2023

Balance: Sartorius -
LE244S-OCE

DE-M