

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 011 005

Volumetric flask, 100 ml, NS 14, A, KB

Batch: 04-20


Serial:

Measuring Nr. 35791

Batch size:

57

Volume:	100,000 ml
Tolerance:	$\pm 0,0500$ ml
Average value:	99,9967 ml
Deviation:	0,0087 ml
Precision (CV):	0,0087 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

31.01.2020

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M