

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 012 BL USP

Vol. flask, 10 ml, NS 10, blue grad., USP, A, conf. cert.


Batch: 36-19

Serial:

Measuring Nr. 32754

Batch size: 450

Volume:	10,0000 ml
Tolerance:	± 0,0200 ml
Average value:	10,0047 ml
Deviation:	0,0053 ml
Precision (CV):	0,0534 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 77

Date: 03.09.2019

Balance: Sartorius -
LE244S-OCE

DE-M